

SUPPORTS AND ORTHOSES

MEDICAL COMPRESSION STOCKINGS

INSOLES AND SHOES

MEASUREMENT TECHNOLOGY

SUPPORTS AND ORTHOSES

## MalleoLoc® and MalleoTrain® S

*The stability program for  
your ankle.*

**NEW**



Rev. 1 - 12/09\_00000000000000

**Bauerfeind AG**  
Triebeser Straße 16  
D-07937 Zeulenroda-Triebes,  
Germany  
Phone +49 (0) 36628-66-40 00  
Fax +49 (0) 36628-66-44 99  
E-mail [info@bauerfeind.com](mailto:info@bauerfeind.com)

**USA**  
Bauerfeind USA, Inc.  
3005 Chastain Meadows Parkway  
Suite 700  
Marietta, GA 30066  
Phone 1 800 423 34 05  
Phone (770) 429 83 30  
Fax (770) 429 84 77  
E-mail [info@bauerfeindusa.com](mailto:info@bauerfeindusa.com)

**United Kingdom**  
Bauerfeind UK  
Phyllis House  
229 Bristol Road  
GB-Birmingham B5 7UB  
Phone +44 (0) 121 446 53 53  
Fax +44 (0) 121 446 54 54  
E-mail [info@bauerfeind.co.uk](mailto:info@bauerfeind.co.uk)



# MalleoLoc®

**NEW**



## INDICATIONS

- Post-operative protection after ligamental suturing/reconstruction
- Conservative treatment of severe ankle joint sprains and ligament ruptures
- Early functional treatment for injuries of the lateral malleolar ligaments (and the bifurcate ligaments)
- Chronic ligament insufficiency

## FEATURES

- Optimum treatment quality in active treatment, with no need for any reshaping.
- The anatomical shape of the orthosis. It fits snugly on the foot and is suitable for wearing with and without a shoe.
- Comfortable, anti-fatigue cushioning as well as a Velcro strap system to wind in the shape of a figure of 8. This guarantees maximum wearing comfort.

# MalleoTrain® S

**NEW**



## INDICATIONS

- Ligament insufficiency
- Slight sprains
- Supination prophylaxis particularly in connection with sports activities
- Post-operative rehabilitation

## FEATURES

- 3-level strap system (elastic, inelastic, semi-elastic) stabilizes the ankle joint when wound in a figure of 8 pattern.
- The securing effect on the ligaments can be individually adjusted depending on the medical indication.