

**Product catalog**  
Medical compression  
stockings

**VenoTrain®**



## Markings on packaging

### Top bands




Silicone dots



Silicone dot pattern



Lace



Honeycomb





### Toe



Closed



Open

					
<i>micro</i>	CREAM				
	HONEY				
	CARAMEL				
	AMARETTO				
	ALMOND				
	ESPRESSO				
	BLACK	●	●	●	●
	GRAPHITE				
	NAVY				
	LAVENDER				
	SILVER				
	DENIM				
RUBY					
<i>micro balance</i>	CREAM				
	CARAMEL	●	●	●	
<i>look</i>	BLACK				
	CARAMEL				
	TOFFEE				
<i>soft</i>	BLACK	●	●	●	●
	GRAPHITE				
	NATURAL				
	CARAMEL				
<i>business</i>	ESPRESSO	●	●	●	● (custom-made)
	BLACK				
	NAVY				
	WHITE				
	TERRA				
<i>sport</i>	BLACK	●			
	TITANIUM				
	WHITE				
<i>impuls+</i>	BLACK	●	●	● (custom-made)	
	CARAMEL				
<i>delight</i>	CARAMEL	a variety of styles			
<i>ulcertec</i>	NATURAL	●			
<i>curaflow</i>	CREAM	a variety of styles			
	CARAMEL	a variety of styles			
	BLACK	a variety of styles			

## Products and health care concepts:

### For better quality of life

As one of the world's leading manufacturers of medical products, Bauerfeind AG has, through its products and service concepts, been making a significant contribution toward maintaining and restoring health and increasing well-being for more than 80 years – for a better quality of life.

Our specialists work with physicians and scientists from all over the world and are fully aware of what is expected of them, in terms of both the tradition and reputation of a family-run business and the high performance levels demanded by today's healthcare market.

Numerous patents and international awards underline our company's flair for innovation and our outstanding expertise in technology and manufacturing. We intend to continue delivering quality "Made in Germany".

#### Supports and orthoses

Whether for prevention, treating symptoms or providing post-operative care, physicians and physiotherapists around the world trust our supports and orthoses, which are available in a wide range of styles for all relevant forms of treatment and set the standard for maximizing therapeutic success.

#### Insoles and shoes

In the orthopedic orthosis segment, we set the standard with traditional products and innovative foot orthoses concepts such as TRIactive and ErgoPad. Comfort and therapeutic shoes and pressure relief orthoses round off a comprehensive range for modern medical foot care.

#### Measurement technology

Medical products will only work if they fit perfectly. With this in mind, Bauerfeind has specifically developed digital measurement systems for medical retailers to improve treatment with medical compression stockings and to enable the manufacture of custom-made foot orthoses.

## VenoTrain medical compression stockings:

### Always the right choice

Whether as a preventive measure, for standard phlebological cases, for advanced vein problems or even for the treatment of venous leg ulcers and lymphedema:

VenoTrain medical compression stockings cover all uses/indications associated with modern compression therapy and are also offer clinically proven effectiveness, unique material properties and a patient-oriented modern design.

The indication-driven treatments cover a wide range of round-knit and flat-knit styles in both custom-made and standard versions. The Perfect-Fit sizing system developed by Bauerfeind on the basis of digital leg measurements offers unique flexibility in ensuring that medical compression stockings from the standard product range fit perfectly.

All VenoTrain products satisfy the requirements of the German Medical Devices Act (Medizinproduktegesetz - MPG) and bear the RAL quality mark. The Oeko-Tex Standard 100 is a quality mark that guarantees particularly high levels of skin tolerance.



Measurement technology:

## Guarantees optimal care for your customers



### **BODYTRONIC 200**

Leg health can now be measured quickly and easily thanks to the BODYTRONIC 200, the compact and portable system for measuring vein function.

With wireless transfer to a PC via Bluetooth, the measurement data can be used to devise preventive measures or obtain additional information for proposed medical treatments.



### **BODYTRONIC 600**

Individual measurements for medical compression stockings need to be highly precise if the proposed treatment is to be effective (quality and test specifications RAL-GZ 387, September 2000). The BODYTRONIC 600 measurement system gives you the security of reliable 3D circumference values and a guaranteed fit and can even link up directly with Bauerfeind's online ordering system.



### **Perfect Fit sizing system**

More than 200,000 legs have already been measured using the Image 3D technology, and these provide the basis for the Perfect-Fit system, which can be used in turn to significantly increase the number of your customers treated with standard sizes of medical compression stockings. This reduces delivery times, simplifies logistics and lowers costs while ensuring items fit as well as possible - all of which benefits everybody.

## Service numbers

### USA


Bauerfeind USA, Inc.  
 3005 Chastain Meadows Parkway  
 Suite 700  
 Marietta, GA 30066  
 Phone 1 800 423 3405  
 Phone (770) 429 8330  
 Fax (770) 429 8477  
 E-mail [info@bauerfeindusa.com](mailto:info@bauerfeindusa.com)

### Canada

Galien  
 1200 Bergar  
 Laval, Quebec, Canada  
 H7L 5A2  
 Office Local (450) 629-9889  
 Office Toll Free 1-877-629-9889  
 Fax Local (450) 629-9895  
 Fax Toll Free 1-877-629-9895  
[www.galiengroup.com](http://www.galiengroup.com)

### United Kingdom

Bauerfeind UK  
 Phyllis House  
 229 Bristol Road  
 GB-Birmingham B5 7UB  
 Phone +44 (0) 121 446 53 53  
 Fax +44 (0) 121 446 54 54  
 E-mail [info@bauerfeind.co.uk](mailto:info@bauerfeind.co.uk)



### Croatia

Bauerfeind d.o.o.  
 Goleška 20  
 HR-10020 Zagreb  
 Phone +385 (0) 1 65 42 855  
 Fax +385 (0) 1 65 42 860  
 E-mail [info@bauerfeind.hr](mailto:info@bauerfeind.hr)

### Singapore

Bauerfeind Singapore Pte Ltd  
 Blk 41 Cambridge Road  
 #01-21 Singapore 210041  
 Phone 065 6396 3497  
 Fax 065 6296 5062  
 E-mail [info@bauerfeind.com.sg](mailto:info@bauerfeind.com.sg)

### Australia

MAS Medical Pty Ltd  
 44 Highbury Road  
 Burwood, VIC 3125  
 Phone 03 9830 6529  
 Fax 03 9830 6517  
 E-mail [info@masmedical.com.au](mailto:info@masmedical.com.au)

# VenoTrain<sup>®</sup> *micro*



## Colors

CREAM

HONEY

CARAMEL

AMARETTO

ALMOND

NAVY\*

ESPRESSO

JAZZ

SWING

GRAPHITE

BLACK

LAVENDER

SILVER

DENIM

RUBY

## Qualities

The internationally popular compression stocking with a particularly high microfiber content for lightweight wearing comfort.

In a wide range of styles, discreet colors and patterns, or more trendy colors to suit any look.

- 50 % microfibers: breathable, kind to the skin and easy to clean
- figure-enhancing panty area (pantyhose) and comfort sole
- opaque full-fashioned stocking quality

### Material content

- Ccl1 65 % polyamide  
35 % elastane
- Ccl2 55 % polyamide  
45 % elastane

### Care instructions

104 °F     

Ccl1 | Ccl2

\* No pattern

Colors may differ due to printing.

## Styles | Standard

	Size	Compression cl. (Ccl)		Edging/Top band				Toe	
		1	2	Top band	Lace	Silicone dots	Silicone dot pattern	Open	Closed
AD	XS (caramel)		●	●					●
	S-XL	●	●	●		●		●*	●
AG	S-XL	●	●	●	●	●	●	●*	●
AT	S-XL	●	●	●				●*	●
ATU	S-XL	●	●	●				●*	●

## Styles | Custom-made solutions

		Compression cl. (Ccl)		Edging/Top band				Toe		
		1	2	Top band	Lace	Honey-comb	Silicone dots	Silicone dot pattern	Open	Closed
AD		●	●	●			●		●	●
AF		●	●	●	●	●	●	●	●	●
AG		●	●	●	●	●	●	●	●	●
AT		●	●	●					●	●
ATU		●	●	●					●	●



# VenoTrain<sup>®</sup> *micro balance*



## Colors

CREAM

CARAMEL

BLACK

## Qualities

Specially designed for sensitive legs, with active skin care substances and reload emulsion:

Simply spray the emulsion onto the stocking for lasting skin care - keeping the skin smooth and supple.

- 50 % microfibers: breathable, kind to the skin and easy to clean
- opaque full-fashioned stocking quality
- stockings are easy to put on and take off

### Material content

- |      |                |
|------|----------------|
| Ccl1 | 65 % polyamide |
|      | 35 % elastane  |
| Ccl2 | 55 % polyamide |
|      | 45 % elastane  |

### Care instructions

104°F

*For treatment with reload emulsion, please see separate product information.*




Ccl1 | Ccl2






*Reload emulsion with active ingredients - lasting skin care, simply by spraying your stockings*

*Colors may differ due to printing.*

## Styles | Standard

	Size	Compression cl. (Ccl)		Edging/Top band				Toe	
		1	2	Top band	Lace	Silicone dots	Silicone dot pattern	Open	Closed
 AD	S-XL	●	●	●				●	●
 AG	S-XL	●	●		●	●	●	●	●
 AT	S-XL	●	●	●				●	●

## Styles | Custom-made solutions

		Compression cl. (Ccl)		Edging/Top band				Toe	
		1	2	Top band	Lace	Silicone dots	Silicone dot pattern	Open	Closed
 AD		●	●	●				●	●
 AG		●	●		●	●	●	●	●
 AT		●	●	●				●	●



# VenoTrain<sup>®</sup> soft | soft S

## Colors

NATURAL


CARAMEL

ESPRESSO

BLACK

NAVY

Ccl 1 | Ccl 2



## Qualities

Soft and robust, with a range of special designs available. In "natural" color version with express manufacture of made-to-measure sizes, ideal for fast post-operative treatment - shipment on the same day for orders before 10 a.m.

For patients with restricted mobility, the VenoTrain soft S (standard design in the "caramel" color version) offers a generous heel area and an especially comfortable instep.

- for him and her
- with a beneficial massage action
- also available in a non-colored ("natural") version

### Material content

Ccl 1 70 % polyamide  
30 % elastane

Ccl 2 65 % polyamide  
35 % elastane

### Care instructions

104 °F

\* VenoTrain soft; see the price list for VenoTrain soft S styles

● Only available with 3 cm silicone dotted band and in Ccl 2.

● Silicone dotted band 3 cm

Colors may differ due to printing.



## Styles | Standard

	Size	Compression cl. (Ccl)		Edging/Top band				Toe	
		1	2	Top band	Honey-comb	Silicone dots	Garter belt	Open	Closed
soft	AD XS		●	●				●	
	AD S-XL	●	●	●		○		●	●
	AG S-XL	●	●	●	●	●	●	●	●
	AT S-XL	●	●	●				●	●
soft S (caramel)	AD XS-XL	●	●			○		●	
	AG S-XL	●	●			●		●	

## Styles | Custom-made solutions\*

	Compr. cl. (Ccl)		Edging/Top band					Toe				Panty-hose top
	1	2	Top band	Lace	Honey-comb	Silicone dots	Garter belt	Open	Closed	No foot	Stump	Open crotch
AD	●	●	●			●		●	●	●	●	
AF	●	●	●	●	●	●		●	●	●	●	
AG	●	●	●	●	●	●	●	●	●	●	●	
AT	●	●	●					●	●	●	●	●
ATU	●	●	●					●	●	●		
ATH	●	●	●					●	●	●		
ATX	●	●	●					●	●	●		
ATE	●	●	●					●	●	●		



# VenoTrain<sup>®</sup> look

## Colors

CARAMEL

TOFFEE

BLACK

GRAPHITE

Ccl1 | Ccl2

## Qualities

A real highlight thanks to its glossy look: figure-enhancing with a push-up effect, stretch comfort and trilobal fiber for a perfect fit.

The fashionable accessory that is also suitable for preventive care and during pregnancy.

- fashionable and elegant glossy effect
- slightly transparent

### Material content





Ccl1	60 % polyamide
	40 % elastane
Ccl2	55 % polyamide
	45 % elastane

### Care instructions





104 °F     

*Colors may differ due to printing.*

## Styles | Standard

	Size	Compression cl. (Ccl)		Edging/Top band		Toe
		1	2	Top band	Silicone dot pattern	
 AD	S-XL	●	●	●		●
 AG	S-XL	●	●		●	●
 AT	S-XL	●	●	●		●
 ATU	S-XL	●	●	●		●

## Styles | Custom-made solutions

	Compression cl. (Ccl)		Edging/Top band			Toe
	1	2	Top band	Lace	Silicone dot pattern	
 AD	●	●	●			●
 AG	●	●	●	●	●	●
 AT	●	●	●			●
 ATU	●	●	●			●



# VenoTrain<sup>®</sup> *discrétion*

## Colors

CLAIR

BEIGE

ESPRESSO

BLACK

JAVA

Classe 1 | Classe 2

## Qualities

The exceptionally fine material of these compression stockings gives them a stylish transparency and lightness.

Available in various designs and versatile colors, these stockings will give you an extra sense of vitality and well-being in your work, leisure and everyday life.

- highly elastic, making the stockings easy to put on
- shapes the legs without constricting them
- breathable fine knit for a long-lasting wearing comfort

### Material content

Classe 1 75 % polyamide  
25 % elastane

Classe 2 70 % polyamide  
30 % elastane

### Care instructions

104 °F     

*VenoTrain discrétion is only available in selected countries.*

*Colors may differ due to printing.*

## Designs | Standard

	Size	Compression		Colors					Toe	
		10-15 mmHg	15-20 mmHg	Clair	Beige	Espresso	Black	Java	Open	Closed
AD Knee stockings	1-4 short et long	●			●	●	●	●	●	●
	1-4 short et long		●	●	●	●	●	●	●	●
AG Thigh-high stockings with silicone lace band	1-4 short et long	●			●	●	●	●	●	●
	1-4 short et long		●	●	●	●	●	●	●	●
AG Thigh-high stockings with silicone dotted band	1-4 short et long	●			●	●	●	●	●	●
	1-4 short et long		●	●	●	●	●	●	●	●
AT Pantyhose	1-4 short et long	●			●	●	●	●	●	●
	1-4 short et long		●		●	●	●	●	●	●

Size	1	2	3	4
cG	42-56	49-63	56-70	63-77
cC	28-36	32-40	35-43	38-46
cB	18-21	21-24	24-27	27-31

Size	short	long
AG   AT	<70	>70
AD	<38	>38

## Designs | Custom-made solutions

	Size	Compression		Colors			Toe	
		10-15 mmHg	15-20 mmHg	Beige	Espresso	Black	Open	Closed
AD Knee stockings	1-4 short et long	●		●	●	●	●	●
	1-4 short et long		●	●	●	●	●	●
AG Thigh-high stockings with silicone lace band	1-4 short et long	●		●	●	●	●	●
	1-4 short et long		●	●	●	●	●	●
AG Thigh-high stockings with silicone dotted band	1-4 short et long	●		●	●	●	●	●
	1-4 short et long		●	●	●	●	●	●
AT Pantyhose	1-4 short et long	●		●	●	●		●
	1-4 short et long		●	●	●	●		●



# VenoTrain<sup>®</sup> act

## Colors

BIANCO

CHIARO

ORO

DAINO

VISONE

CAFFÉ

NERO

FUMO

BLU

70 | 140 | 280 den

## Qualities

High-quality support stockings for men and women, designed to provide beneficial relief for the legs and prevent vein problems.

Available in various designs and colors, in opaque or transparent versions and with three compression strengths.

- wide range of attractive products
- six standard sizes for an optimum fit

### Material content

70 den AD 82 % polyamide  
18 % elastane  
AG 81 % polyamide  
19 % elastane  
AT 78 % polyamide  
22 % elastane

140 den AD 62 % polyamide  
38 % elastane  
AG 72 % polyamide  
28 % elastane  
AT 74 % polyamide  
26 % elastane

280 den AD 37 % polyamide  
11 % elastane  
52 % cotton

### Care instructions

86°F     

*VenoTrain act is only available in selected countries.*

*Colors may differ due to printing.*



## Designs | Standard

			Colors								
			Bianco	Chiario	Oro	Daino	Visone	Caffé	Nero	Fumo	Blu
	70 den smooth	1-5		●	●	●	●	●	●		
	140 den mesh	1-6		●	●	●	●	●	●		
	70 den smooth	2-5		●	●	●	●	●	●		
	140 den mesh	2-6		●	●	●	●	●	●		
	140 den mesh	1-5			●				●		
	70 den patterned	1-5							●		
	70 den smooth	1-5		●	●	●	●	●	●		
	140 den mesh	1-6		●	●	●	●	●	●		
	70 den smooth	1-5		●	●	●	●	●	●		
	140 den mesh	1-5		●	●	●	●	●	●		
	280 den	1-3	●					●	●	●	●

VenoTrain® act

		Pantyhose   Thigh-high stockings (AG   AT   ATU)												
		Weight	40	45	50	55	60	65	70	75	80	85	90	+ kg
Height	145 cm													
	150 cm	1												
	155 cm			2										
	160 cm				3									
	165 cm					4								
	170 cm						5							
	175 cm							6						
	185 cm													
			99	110	121	132	143	154	165	176	187	198	209	+ lbs

Knee stockings (AD)					
Size	1	2	3	4	5
Shoe size	35-37	38-39	40-41	42-43	>43

Knee stockings (AD) - cotton			
Size	1	2	3
Shoe size	39-40	41-42	>42



# VenoTrain<sup>®</sup> *business*

## Colors

WHITE

TERRA

BLACK

GRAPHITE

NAVY

Ccl 1 | Ccl 2



## Qualities

The men's stocking made of cotton and Tactel<sup>®</sup> sports fiber offers first-class wearing comfort with thermal regulation.

Highly elastic LYCRA<sup>®</sup> soft and the extra-wide supportive band ensure the item is easy to put on and stays perfectly in position.

- high-quality men's stocking with a classically discreet ribbed structure
- innovative mix of materials
- reinforced sole
- Timbrelle<sup>®</sup> high-tech fiber for long-lasting freshness around the foot

### Material content


Ccl 1	80 % polyamide (incl. 50 % Tactel)
	14 % elastane
	6 % cotton
Ccl 2	71 % polyamide (incl. 45 % Tactel)
	24 % elastane
	5 % cotton

### Care instructions


104 °F     

Colors may differ due to printing.

## Styles | Standard

	Size	Compression cl. (Ccl)		Edging/Top band	Toe
		1	2		
 AD	S-XL	●	●	●	●

## Styles | Custom-made solutions

	Compression cl. (Ccl)	Edging/Top band	Toe
 AD	●	●	●

# VenoTrain<sup>®</sup> sport



## Colors

WHITE

TITANIUM

BLACK

20 mmHg



## Qualities

Specifically tailored to athletes' requirements for improving endurance and promoting regeneration.

Muscle-active zones in the knitted fabric stabilize and support the calf muscle.

Optimal temperature regulation thanks to the very high microfiber content.

- 50 % microfibers: breathable, kind to the skin and easy to clean
- reinforced antibacterial sole
- outstanding medical effectiveness
- four standard sizes, each in two lengths

### Material content


65 % polyamide  
35 % elastane

### Care instructions

104 °F     

Colors may differ due to printing.

## Styles | Standard

	<i>Size</i>	<i>Compression</i>	<i>Edging/Top band</i>	<i>Toe</i>
			Top band	Closed
 AD	S-XL	20 mmHg	●	●

# VenoTrain® *impuls*<sup>+</sup>



## Colors

CARAMEL

BLACK

Ccl2 | Ccl3



## Qualities

The high working pressure, even during small muscle movements, produces a really deep action. Easy to put on and breathable thanks to microfibers. Express manufacture of made-to-measure sizes.

Highly effective for advanced vein problems, post-thrombotic syndrome, providing follow-up care for healed ulcers and providing therapy in the early stages of lymphedema.

- short-stretch material with skin-friendly microfiber
- highly differentiated V6+ sizing system

### Material content

Ccl2 80 % polyamide - including 44 % microfiber  
20 % elastane



Ccl3 75 % polyamide - including 42 % microfiber  
25 % elastane

### Care instructions






86°F

● Silicone dotted band 3 cm  
Colors may differ due to printing.

## Styles | Standard

	Size	Compression cl. (Ccl)		Edging/Top band			Toe	
		2	3	Top band	Silicone dots	Garter belt	Open	Closed
 AD	I-VI	●	●	●	◐		●	●
 AG	I-VI	●	●	●	●	●	●	●

## Styles | Custom-made solutions

		Compression cl. (Ccl)		Edging/Top band			Toe	
		2	3	Top band	Silicone dots	Garter belt	Open	Closed
 AD		●	●	●	●		●	●
 AG		●	●	●	●	●	●	●
 AT		●	●	●			●	●
 ATH		●	●	●			●	●
 ATX		●	●	●			●	●



# VenoTrain<sup>®</sup> *ulcertec*



## Colors

NATURAL

Overstocking

WHITE

Liner

34-39 mmHg  
40-46 mmHg



## Qualities

The combination of an overstocking with a patented rhomboid knit and a skin-friendly liner for fixing the wound dressing in place makes the item easy to put on without any help from others.

The liners are also designed to provide a continuous low level of compression during the night, which can be individually adjusted for the set of three (at 10 or 18 mmHg), and are suitable for washing at high temperatures.

- easy to put on
- constant pressure
- discreet appearance
- X5 sizing system

### Material content of overstocking

77% polyamide  
23% elastane

### Material content of liner

10 mmHg 89% polyamide  
11% elastane  
18 mmHg 72% polyamide  
28% elastane

### Care instructions for overstocking

104°F     

### Care instructions for liner

203°F     

Colors may differ due to printing.

## Styles | Standard

Size		Compression	Edging/Top band	Toe	
			Top band	Open	Closed
 AD	XS-XL	34-39 mmHg and 40-46 mmHg	●	Overstocking	Liner

### VenoTrain ulcertec 39 Liner 10 mmHg Overstocking 23-27 mmHg

- compression strength 34-39 mmHg
- non-problematic venous ulcer of the leg with good tendency to heal
- light liner for patients who can only tolerate a little compression at night

### VenoTrain ulcertec 46 Liner 18 mmHg Overstocking 23-27 mmHg

- compression strength 40-46 mmHg
- treatment-resistant ulcers
- ulcers with a poor tendency to heal
- with tolerance and a good level of acceptance of compression

# VenoTrain® *delight*



## Colors

CARAMEL

Ccl1 - Ccl4



## Qualities

Available as a flat-knit arm sleeve or leg stocking with skin-friendly, latex-free materials and exceptional fit. The special gripRelief® structure at the top ensures the item holds comfortably and securely in place.

For the treatment of advanced venous and lymphatic symptoms

- breathable double-stretch knit
- extra-flat seams
- wide variety of styles

### Material content - stockings

Ccl2	80 % polyamide 20 % elastane
Ccl3	79 % polyamide 21 % elastane
Ccl4	75 % polyamide 25 % elastane

### Material content - arm sleeves

Ccl1 2	80 % polyamide 20 % elastane
--------	---------------------------------

### Care instructions

86°F

● *Optionally with T-heel or pads*

*Please refer to the current price list for further styles.*

*Colors may differ due to printing.*

## Styles | Custom-made solutions

	Compression class (Ccl)				Edging/Top band			Toe			Additional features	
	1	2	3	4	Top band	Silicone dots	Garter belt	Open	Closed	Diagonal cut	T-heel	Visco-elastic pads
AD		●	●	●	●	●		●	●	●	◐	◐
AF		●	●	●	●	●		●	●	●	◐	◐
AG		●	●	●	●	●	●	●	●	●	◐	◐
AT		●	●		●			●	●	●	●	
ET		●	●		●							
CT		●	●		●							
CG	●	●			●	●						
CG1	●	●			●	●						

# VenoTrain<sup>®</sup> curaflow



## Colors

CREAM

CARAMEL

BLACK



Ccl1 - Ccl4

## Qualities

The high microfiber content and the innovative way these are processed make this item unique. Extra-flat sensitive stitching and special comfort zones ensure this flat knit treatment provides excellent wearing comfort.

Available for the leg, foot, arm, hand and thorax.

- smooth microfibers make the item easy to put on
- breathable knit
- extra-soft behind the knee for problem-free mobility
- special shapes to suit individual requirements

### Material content

Ccl1-4 75% polyamide  
25% elastane

### Care instructions

86°F     

Please refer to the VenoTrain curaflow product catalog or the price list for additional styles (hand, foot and toe).

Colors may differ due to printing.

## Styles

### Custom-made solutions, leg

	Compression class (Ccl)			
	1	2	3	4
AA1	●	●	●	
AD	●	●	●	●
AG	●	●	●	●
AT	●	●	●	●
ET	●	●	●	●
FT	●	●	●	●
CT	●	●	●	●

### Custom-made solutions, arm

	Compression cl. (Ccl)			Edging		
	1	2	3	Straight	Sloping	Without hand
CG	●	●	●		●	●
AG	●	●	●		●	
AF	●	●	●	●		
AE	●	●	●	●		
AD	●	●	●			
AC1	●	●	●			

### Custom-made solutions, thorax

	Compression class (Ccl)
	1
Body	●
T-shirt	●
Bolero under the bust	●

## Standard product range | round-knit products

Specifications		Size	Compression class (Ccl)			Edging/Top band						Toe	
			1	2	3	Top band	Lace	Honeycomb	Silicone dots	Silicone dot pattern	Garter belt	Open	Closed
<b>VenoTrain</b>													
<b>micro</b>	<b>AD</b>	XS		●		●							●
		S-XL	●	●		●			●			●	●
	<b>AG</b>	S-XL	●	●		●*	●		●	●		●	●
	<b>AT</b>	S-XL	●	●		●						●	●
	<b>ATU</b>	S-XL	●	●		●						●	●
<b>micro balance</b>	<b>AD</b>	S-XL	●	●		●						●	●
	<b>AG</b>	S-XL	●	●			●		●	●		●	●
	<b>AT</b>	S-XL	●	●		●						●	●
<b>soft</b>	<b>AD</b>	XS		●		●						●	
		S-XL	●	●		●			○			●	●
	<b>AG</b>	S-XL	●	●		●*	●			●*		●	●
	<b>AT*</b>	S-XL	●	●		●						●	●
<b>soft S</b>	<b>AD</b>	XS-XL	●	●					○			●	
	<b>AG</b>	XS-XL	●	●					●			●	
<b>look</b>	<b>AD</b>	S-XL	●	●		●							●
	<b>AG</b>	S-XL	●	●						●			●
	<b>AT</b>	S-XL	●	●		●							●
	<b>ATU</b>	S-XL	●	●		●							●
<b>business</b>	<b>AD</b>	S-XL	●	●		●							●
<b>sport</b>	<b>AD</b>	S-XL	20 mmHg			●							●
<b>impuls*</b>	<b>AD</b>	S-XL		●	●	●			○			●	●
	<b>AG</b>	S-XL		●	●	●			●		●	●	●
<b>ulcertec</b>	<b>AD</b>	XS-XL	34-39 mmHg and 40-46 mmHg			●						Over- stocking	Liner

XS sizes and VenoTrain soft S are only available in color caramel

○ Only available with 3 cm silicone dotted band and in Ccl 2.

○ 3 cm silicone dotted band

## Custom product range | round-knit products

Specifications		Compression cl. (Ccl)			Abschluss/Haftband						Edging/Top band				Pantyhose top
		1	2	3	Top band	Lace	Honeycomb	Silicone dots	Silicone dot pattern	Garter belt	Open	Closed	No foot	Stump	Open crotch
VenoTrain															
micro	AD	•	•		•			•			•	•			
	AF	•	•		•	•	•	•	•		•	•			
	AG	•	•		•	•	•	•	•		•	•			
	AT	•	•		•						•	•			
	ATU	•	•		•						•	•			
micro balance	AD*	•	•		•						•	•			
	AG	•	•			•		•	•		•	•			
	AT	•	•		•						•	•			
soft	AD	•	•		•			•			•	•	•	•	
	AF	•	•		•	•	•	•			•	•	•	•	
	AG	•	•		•	•	•	•		•	•	•	•	•	
	AT	•	•		•						•	•	•	•	•
	ATU	•	•		•						•	•			
	ATH	•	•		•						•	•			
	ATX	•	•		•						•	•			
	ATE	•	•		•						•	•			
soft S	AD	•	•		•			•			•	•			
	AG	•	•				•	•			•	•			
look	AD	•	•		•							•			
	AG	•	•		•	•			•			•			
	AT	•	•		•							•			
	ATU	•	•		•							•			
business	AD	•	•		•						•				
impuls*	AD		•	•	•			•			•	•			
	AG		•	•	•			•		•	•	•			
	AT		•	•	•						•	•			
	ATH		•	•	•						•	•			
	ATX		•	•	•						•	•			

Custom product range | flat-knit products – please refer to VenoTrain delight/curaflow





## Sizing system Perfect-Fit | round-knit products

VenoTrain look, VenoTrain micro, VenoTrain micro balance, VenoTrain soft, VenoTrain soft S, VenoTrain business

Circumference	XS	S		M		L		XL	
	Normal	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus
cG top band, wide*		51-57	59-66	56-63	62-69	60-67	66-74	63-71	71-80
cG top band, standard	44-50	45-51	52-59	49-56	55-62	53-60	58-66	55-63	62-71
cG		45-57	52-66	49-63	55-69	53-67	58-74	55-71	62-80
cF	39-50	40-52	43-55	43-55	46-60	45-59	50-64	49-63	52-68
cE	32-38	33-40	36-44	35-43	39-49	38-48	42-53	42-53	47-58
cD	28-31	29-35	32-38	32-38	35-42	34-42	37-45	37-45	41-49
cC	29-33	31-37	34-40	34-40	37-44	37-45	39-47	40-48	43-51
cB1	23-26	25-29	27-31	27-32	29-34	29-35	31-37	32-38	34-40
cB	18-20	20-22		22-25		25-28		28-31	
cY	27-31	29-33		31-35		34-40		37-44	
cA	17-22	19-24		22-27		25-31		28-35	

Length	short	long	Type
ℓD	35-40	40-46	AD
ℓG	61-70	70-80	AG, AT

\* All products except of VenoTrain soft S and VenoTrain business


## Sizing system | VenoTrain sport

Circumference	S	M	L	XL
cC	31-40	34-44	37-47	40-51
cB	20-22	22-25	25-28	28-31

Length	short	long	Type
ℓD	35-40	40-46	AD

Indication of size in cm


## Sizing system V6+ | VenoTrain impuls+



Circumference	I		II		III		IV		V		VI	
	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus	Normal	Plus
cG	45-53	51-59	48-56	54-62	51-59	57-66	54-63	60-69	57-66	63-73	60-70	66-76
cF	39-45	43-50	42-49	46-53	45-52	49-57	48-56	52-60	51-59	55-64	54-63	58-67
cE	32-37	35-40	34-39	37-42	36-41	39-44	38-44	41-47	40-46	43-49	42-48	45-51
cD	28-32	30-34	30-34	32-37	32-37	34-39	34-39	36-41	36-41	38-43	38-43	40-46
cC	29-33	31-35	31-35	33-38	33-38	35-40	35-40	37-42	37-42	39-44	39-44	41-47
cB1	23-27	24-28	25-29	26-30	27-31	28-32	29-34	30-35	31-36	32-37	33-38	34-39
cB	19-20,5		20,5-22		22-24		24-26		26-28,5		28,5-31	
cY	25-29		26-31		28-33		30-35		32-38		35-41	
cA	19-22		20-23		22-25		24-27		26-30		28-33	

Length	short	long	Type
ℓD	35-39	40-44	AD
ℓG	62-70	71-80	AG

## Sizing system X<sup>5</sup> | VenoTrain ulcertec



Circumference	XS	S	M	L	XL
cC	28-35	31-39	34-42	37-47	41-52
cB	18-20	20-23	23-26	26-30	30-34

Length	short	long	Type
ℓD	bis 40	ab 40	AD

# VenoTrain® glider

VenoTrain glider is the practical aid for putting on and taking off all types of medical compression stockings.

VenoTrain glider can be prescribed by your physician.

## Hard-wearing material

- stable and tear-resistant, extremely lightweight

## Saves time and space

- can be put on and taken off quickly and easily
- also ideal when traveling



# VenoTrain® Accessories



A **lingerie bag** is a safe way to protect VenoTrain medical compression stockings in every machine wash.



**Rubber gloves** make donning compression stockings a snap.



The **special detergent** cleans gently, ensuring careful treatment and providing lasting protection for your products.



The **non-skid sole** prevents slipping out.



The **silicone top band cleaning cloth** ensures the top band keeps the stockings in place, even during the longest days.




Doing simple exercises with the **exercise band** after sitting or standing for long periods of time revitalizes your legs in no time.


# A systematic approach to medical compression stockings: For treatment suited to the individual indication

Diagnostic findings	No visible signs of a venous disorder	Spider veins Reticular varices	Varicosis Imminent functional insufficiency	Varicosis with edema Functional insufficiency	Varicosis with trophic skin changes Atrophy blanche, dermatoliposclerosis	Healed ulcus cruris venosum	Active ulcus cruris venosum
CEAP-classification*	<b>C0</b>	<b>C1</b>	<b>C2</b>	<b>C3</b>	<b>C4</b>	<b>C5</b>	<b>C6</b>

Medical condition	Mild discomfort <b>phlebological standard case</b>	Severe discomfort <b>phlebological severe case</b>
-------------------	---	---



Compression classes	<b>Ccl 1</b> 18-21 mmHg	<b>Ccl 2</b> 23-32 mmHg	<b>Ccl 3</b> 34-46 mmHg
---------------------	----------------------------	----------------------------	----------------------------



\* Partsch, H.: Klassifikation von chronischen Venenerkrankungen heute, Vasomed 1/2003  
\*\* E. Rabe et al. Bonner Venenstudie der DGP, Phlebologie 1/2003

MEDICAL COMPRESSION STOCKINGS

SUPPORTS AND ORTHOSES

INSOLES AND SHOES

MEASUREMENT TECHNOLOGY

#### **Australia**

MAS Medical Pty Ltd  
44 Highbury Road  
Burwood, VIC 3125  
Phone 03 9830 6529  
Fax 03 9830 6517  
E-mail [info@masmedical.com.au](mailto:info@masmedical.com.au)

#### **Canada**

Galien  
1200 Bergar  
Laval, Quebec, Canada  
H7L 5A2  
Office Local (450) 629-9889  
Office Toll Free 1-877-629-9889  
Fax Local (450) 629-9895  
Fax Toll Free 1-877-629-9895  
[www.galiengroup.com](http://www.galiengroup.com)

#### **Singapore**

Bauerfeind Singapore Pte Ltd  
Blk 41 Cambridge Road  
#01-21 Singapore 210041  
Phone 065 6396 3497  
Fax 065 6296 5062  
E-mail [info@bauerfeind.com.sg](mailto:info@bauerfeind.com.sg)

#### **United Kingdom**

Bauerfeind UK  
Phyllis House  
229 Bristol Road  
GB-Birmingham B5 7UB  
Phone +44 (0) 121 446 53 53  
Fax +44 (0) 121 446 54 54  
E-mail [info@bauerfeind.co.uk](mailto:info@bauerfeind.co.uk)

#### **USA**

Bauerfeind USA, Inc.  
3005 Chastain Meadows Parkway  
Suite 700  
Marietta, GA 30066  
Phone 1800 423 3405  
Phone (770) 429 8330  
Fax (770) 429 8477  
E-mail [info@bauerfeindusa.com](mailto:info@bauerfeindusa.com)

Rev. 9 – 01/13\_27120101001001

#### **Bauerfeind AG**

Triebeser Straße 16  
07937 Zeulenroda-Triebes  
Phone +49 (0) 36628-66-2000  
Fax +49 (0) 36628-66-2999  
E-mail [info@bauerfeind.com](mailto:info@bauerfeind.com)